

# Pulseq Principles

Maxim Zaitsev

*Division of Medical Physics, Dept. of Radiology,  
University Medical Center Freiburg, Germany*

*Feb 8, 2024*

# What is *Pulseq*?

---

- *Pulseq* is a language to describe MR pulse sequences
- *Pulseq* sequences are fixed successions of RF and gradient pulses and ADC events



- *Pulseq* is the software to generate such pulse sequence descriptions
- *Pulseq* scripts can re-generate *Pulseq* sequences to accommodate user input

*Pulseq* ecosystem includes sequences, software to generate them and software and hardware to consume them


# *Pulseq* Philosophy

---

- Minimize effort for implementation and support on hardware
  - Lean sequence-to-hardware interface
- Remove the thresholds in sequence programming
  - Make simple things truly simple
- Make researcher-oriented features accessible
  - Arbitrary gradients, arbitrary RF, free ordering, X-nuclei, ...
- Prevent typical sources of (human) errors
  - Avoid timing errors with “overlapping” gradients
  - Make data flag and counter setting optional/unnecessary
- **Promote open-source thinking, sharing and exchange!**



# Pulse sequence definition in *Pulseseq*



- Block 1:  
gradient and RF
- Block 2:  
only gradients
- Block 3:  
gradient and ADC
- Block 4:  
only gradients
- Block 5:  
gradient and RF ...

- Sequence is a concatenation of non-overlapping blocks
- Gradients do not have to start or end at 0 at the block boundaries

# *Pulseq* block concept in detail

---

- Each block may contain following events:
  - One optional gradient pulse per axis
  - One optional RF pulse
  - One optional ADC event
- Individual events may define own start delays
- All events in the block overlap in time
- Duration of the block is defined by the longest event
  - Matlab/Python toolboxes use “dummy” delay objects to make blocks longer
- Explicit sequence description
  - No loops, no dependent parameters – like a recorded piece of music!

# High-level programming environments

---

- Matlab *Pulseq* toolbox
- Python *PyPulseq* toolbox



- Further options
  - TOPPE is primarily targeted at GE but can import and export *pulseq* files (Jon-Fredrik Nielsen will talk about it today)
  - GammaStar can export *pulseq* files
  - JEMRIS Bloch simulator can export *pulseq* files
  - CoreMRI Bloch simulator can export *pulseq* files
  - ...

# Matlab *Pulseseq* workflow

```
system = mr.opts('MaxGrad',30,'GradUnit','mT/m',...
    'MaxSlew',170,'SlewUnit','T/m/s');
seq=mr.Sequence(system);

fov = 220e-3; Nx=64; Ny=64; TE = 10e-3; TR = 20e-3;

[rf, gz] = mr.makeSincPulse(15*pi/180,system,'Duration',4e-3,...
    'SliceThickness',5e-3,'apodization',0.5,'timeBwProduct',4);

gx = mr.makeTrapezoid('x',system,'FlatArea',Nx/fov,'FlatTime',6.4e-3);
adc = mr.makeAdc(Nx,'Duration',gx.flatTime,'Delay',gx.riseTime);
gxPre = mr.makeTrapezoid('x',system,'Area',-gx.area/2,'Duration',2e-3);
gzReph = mr.makeTrapezoid('z',system,'Area',-gz.area/2,'Duration',2e-3);
phaseAreas = ((0:Ny-1)-Ny/2)*1/fov;

delayTE = TE - mr.calcDuration(gxPre) - mr.calcDuration(rf)/2 ...
    - mr.calcDuration(gx)/2;
delayTR = TR - mr.calcDuration(gxPre) - mr.calcDuration(rf) ...
    - mr.calcDuration(gx) - delayTE;
delay1 = mr.makeDelay(delayTE);
delay2 = mr.makeDelay(delayTR);

for i=1:Ny
    seq.addBlock(rf,gz);
    gyPre = mr.makeTrapezoid('y',system,'Area',phaseAreas(i),...
        'Duration',2e-3);
    seq.addBlock(gxPre,gyPre,gzReph);
    seq.addBlock(delay1);
    seq.addBlock(gx,adc);
    seq.addBlock(delay2)
end

seq.write('gre.seq')
```

*a runnable gradient echo sequence code  
(similar to Siemens' example miniFlash)*


- Define the system properties
- Define high-level parameters (convenience)
- Define pulses and ADC objects used in the sequence
- Calculate the delays and reordering tables
- Loop and define sequence blocks
- Duration of each block is defined by the duration of the longest event
- *Copy '\*.seq' to the scanner and run it!*



# PyPulseseq workflow

- **PyPulseseq** is a close replica of the original Pulseseq toolbox that does not require a MATLAB license
- Runs in many Python environments, e.g. as notebook in Jupyter (<http://jupyter.org/>) or Google Colaboratory
- Identical workflow:
  - Define the system properties
  - Define high-level parameters (convenience)
  - Define pulses and ADC objects used in the sequence
  - Calculate the delays and reordering tables
  - Loop and define sequence blocks
  - *Download '\*.seq' to the scanner and run it!*

No local installation needed – ideal for education!





# How to design a sequence in Pulseseq

---

## **conceptual design steps**

- Step 1: split the time axis into blocks
- Step 2: assign events to the blocks

## **practical implementation steps**

- Step 3: create/calculate all events
- Step 4: populate the blocks and add them to the sequence

## **validation steps**


- Step 5: check timing, verify k-space trajectory, check hardware and PNS limits, mechanical resonances, etc...

# *Pulseq* objects & blocks

---


- *Pulseq* objects are created and then added to *Pulseq* blocks
  - Gradient & RF pulses  
    mr.makeTrapezoid(...) or pp.make\_trapezoid(...)  
    mr.makeSincPulse(...) or pp.make\_sinc\_pulse(...)
  - ADC objects  
    mr.makeADC(...) or pp.make\_adc(...)
  - Delays, extension objects, etc.  
    (data labels, cardiac trigger directives, trigger pulses)
  - Use seq.addBlock(...) or seq.add\_block(...) to add objects & blocks
- Block concept in *Pulseq* is probably the most demanding part
  - Blocks and objects need to be aligned to raster times

# Example 1: simple gradient echo




- No overlapping gradient ramps on different axes
- Events are clearly separated

# Example 1: basic sequence display options




- Pulseseq 6-panel plot: `seq.plot()`
  - ADC, RF magnitude, RF phase, Gx, Gy, Gz
  - Each event plotted in its own color

# Example 1: display block structure





- Advantageous to separate PE & PR gradients into different blocks
- To visualize block structure:
  - `seq.plot('showBlocks',true,'timeDisp','us');`

# Plot gradient waveforms



- Plot entire waveforms for all axes
  - Native gradient unit in *Pulseseq*: Hz/m


# Basic sequence display options: k-space



- Plot k-space time evolution or 2D (or even 3D) trajectories
  - Native k-space unit in *Pulseq*:  $\text{m}^{-1}$




# Example 2: echo planar imaging (EPI)



- Two RF pulses need to be in separate blocks
- Each ADC event needs to be in a separate block
  - Challenge with blips, readout ramp and optimal sampling window

# EPI block structure




- One possible solution:
  - Keep readout gradient as trapezoid
  - Convert blips to shapes and split them at the center

# *Pulseq* blocks summary

---

- Block concept takes some time to get used to
- Blocks help to organize events and eliminate timing errors
- There is a lot of flexibility
  - Different strategies possible
- Some interpreters expose additional limitations
  - Explicit and implicit delays, number of ADCs per TR, etc...
- **Blocks make it easier for the interpreter to play things out**

# Pulseseq framework overview



# Before you start on the scanner

---

- Peripheral nerve analysis
- Avoiding acoustic resonances (for fast sequences)

# *Pulseq* PNS prediction

---

- Peripheral neural stimulation (PNS) prediction on Siemens is based on the so-called SAFE model

F.X. Herbank and M. Gebhardt. SAFE-Model - A New Method for Predicting Peripheral Nerve Stimulations in MRI. ISMRM 2000, #2007.

<https://cds.ismrm.org/ismrm-2000/PDF7/2007.PDF>

- Open-source implementation by Filip Szczepankiewicz and Thomas Witzel:

[https://github.com/filip-szczepankiewicz/safe\\_pns\\_prediction](https://github.com/filip-szczepankiewicz/safe_pns_prediction)



- Direct interface in Matlab-Pulseq

# Caveat & Solution


---

- Need scanner-specific model parameters
- Extract your gradient system description file from IDEA
  - Type “sys” in the IDEA shell selecting your system  
(you may need to select something else first to see the verbose output)
  - Note the strings after “GPA Type” and “GC Type”
  - Go to `C:\MIDEA\N4_VE####\n4\pkg\MrServers\MrMeasSrv\Config\InitMeas` and pick the file named `MP_GPA_<your_GPA>_<your_coil>.asc` and copy it somewhere where your Matlab or Python can read it
  - *Path under NumarsX is left as an exercise for the interested reader*
- You can now predict PNS in the same way IDEA and scanner do it




# Example: DW-EPI

- $b\text{Factor}=1000$ ,  $TE=78\text{ms}$ ,  $G_{\text{max}}=38\text{mT/m}$   $SR=180\text{T/m/s}$



# PNS Stimulation: DW-EPI

```
seq.calcPNS('idea/asc/MP_GPA_K2309_2250V_951A_AS82.asc');
```




# Acoustic resonance analysis

---

- Extract your gradient system description file from IDEA as above
- Create your sequence in memory  
(populate the “seq” object with events)
- Run “gradSpectrum.m” script in “matlab/demoUnsorted”
  - Set ascName=‘....’ to your gradient system .asc file to see resonances marked

# Acoustic analysis: DW-EPI


---



`seq.sound();` % another option for acoustic assessment in Matlab ;-)

# Pulseq Siemens interpreter

- Just a “normal” sequence
  - Loads its “content” from a Pulseq file
  - Almost all aspects of the sequence are pre-defined in the Pulseq file
  - FOV positioning and rescaling possible
- Based on miniFlash
  - No product code
  - No hacks, no backdoors
- Distributed as a C2P package in **source form**
- Standard SAR calculation
- “Gradient Health” libBalance applicable to all sequences
- PNS and acoustic resonance analysis possible in Matlab
- **Safety equal or higher than a typical IDEA sequence**



# *Pulseq* on Siemens platforms


---

- Over 100 C2P sites
- *Pulseq* works well on:
  - All Numaris4 platforms (tested on vb15...ve12u) and numerous hardware platforms (Symphony, Trio, 7T, 3T Connectom, Skyra, Prisma,...)
  - All released NumarisX versions (xa11, xa20, xa30, xa40, xa5x, xa6x ...)



# Pulseq on Siemens scanners

- *Optional initial step:*  
connect your  
PC to the scanner
- Save the .seq file on the  
scanner as external.seq
- Run the  
interpreter\_sequence  
on the scanner
- *Optional step:*  
stream raw data  
to your PC with  
NIH\_DataCatcher
- or export raw data  
manually







# IDEA simulation with *Pulseseq*



Pulseseq interpreter sequence can also be used with the Siemens' IDEA

1. Save your .seq file as %CustomerSeq%/Pulseseq/external.seq
2. In the IDEA command run **sim**



# *Pulseq* interpreter parameters

- Yes, you can still change few things on the console



# *Pulseq* for Siemens: image recon


---

- Image reconstruction is up to the researcher
- Integration with ICE is possible for some sequences
  - 2D GRE
  - 2D EPI with ramp sampling
  - 3D MPRAGE with GRAPPA
- Online & offline reconstruction with Gadgetron for some sequences
- Offline reconstruction in Matlab & Python
  - Examples for 2D / 3D Cartesian reconstruction
  - Simple gridding reconstruction
  - Example of automated BART reconstruction
- All-in-one FIRE solution : Marten Veldmann et al, MRM 2022, doi:10.1002/mrm.29384

Qingping Chen @  
ISMRM

# Cross-platform sequences with *Pulseq*

- MR physics-oriented workflow
  - Write your sequences from scratch
  - Non-Cartesian readouts, user-defined gradient shapes and custom RF pulses
  - Advanced visualization and analysis tools
  - Automatic k-space calculation
- Pulseq files play out on many scanners
  - Siemens & GE : works
  - Philips : two working interpreters @ISMRM24
  - *Bruker* ???
- New release v1.4.2 (Dec 2023)



## Acknowledgements:

Berkin Bilgic

Frank Zijlstra

Jon-Fredrik Nielsen

Moritz Zaiss

Qiang Liu

Sebastian Littin

Borjan Gagoski

Imam Shaik

Juergen Hennig

Naveen Murthy

Qingping Chen

Will Grissom

Douglas Noll

Jeff Fessler

Mojtaba Shafiekhani

Niklas Wehkamp

Scott Peltier

Yogesh Rathi

# THANK YOU FOR YOUR ATTENTION!



Supported by NIH U24-NS120056 (Nielsen, Zaitsev)  
and R01-EB032378 (Rathi, Bilgic)

