

Echo Planar Imaging (EPI)

Maxim Zaitsev, Ph.D.

University of Freiburg – Medical Centre

Dept. of Radiology – Medical Physics

February 2024

EPI: a historical perspective

- One of the first Fourier imaging techniques^[1]
 - Has immediately polarized the MR community into lovers and haters
- Gained importance only in late 90th of the past century
 - Accessibility of actively-shielded gradients
 - New applications: fMRI, DTI

366

MANSFIELD AND PYKETT





FIG. 9. Pulse and field gradient timing diagram for a two-dimensional echo-planar imaging experiment.


EPI has gone a long way...

EPI of an oil-filled annulus [1]




Sir Peter Mansfield

1.4mm³ EPI at 7T [2]



Maxim Zaitsev

0.45mm³ EPI at 7T Impulse [3]






David Feinberg

[1] Mansfield and Pykett "Biological and medical imaging by NMR" JMR 29,355-73 (1978)

[2] Speck O, Stadler J, Zaitsev M. "High resolution single-shot EPI at 7T" MAGMA 21:73-86 (2008)

[3] Feinberg D et al. Nature Methods 20, 2048-57 (2023)

Today: EPI is a working horse...





A plethora of EPI variants

- Segmented, read-out segmented, PROPELLER-EPI,...
- Accelerated: parallel imaging, multi-band, wave CAIPI,...
- Hybrid gradient/RF encoding (GRASE, RASOR, ...)
- Contrast preparation modules (DTI, etc.)
- Focus on **single-shot EPI**: the most common collection of artifacts


EPI Sequence

- Basic recipe:
 - Oscillating read-out gradient
 - Phase encoding “blips”




Why EPI is special?

- Single shot \neq instantaneous
 - Signal evolution during readout
- Alternating readout polarity
 - Gradient & receiver imperfections
- Fourier reconstruction
 - Perfect encoding?
 - Deviations lead to artifacts...
 - Smooth deviations \rightarrow blurring, distortions
 - Periodic modulation \rightarrow ghosting



EPI artifacts at glance




EPI optimization goals

- Avoid/compensate signal modulation due to oscillating readout
- Maximize sampling efficiency (ADC-on time)
 - Improve signal-to-noise ratio
- Shorten the echo train
 - Reduce blurring due to signal decay
 - Reduce distortions
- Avoid or suppress off-resonance effects
 - Suppress fat
 - Make sure B_0 is as homogeneous as possible
- Use B_0 distortion compensation techniques if possible


Basic EPI Readout in Pulseseq

- Standard 2D excitation
- Define RO duration, calculate ADC and gradient objects
- Calculate PE and RO pre-phasing gradients, PE blips




EPI trajectory video

powered by the pypulseq
Trajectory Animation Tool
by Frank Zijlstra




Get Familiar with EPI Data

- Delays between the gradient and the RF subsystems (RF+ADC)
 - Echoes are shifted in time → in opposite directions in k-space
 - N/2 ghosting artifacts



EPI N/2 Ghosts – 3-echo Reference





- Modulation in K-space \rightarrow “Ghosts” in the image
- Full or partial calibration, image-based methods, coil arrays, data driven...



F. Hennel,
Magma 1999 9:134.
M. Buonocore,
MRM 2001, 45:96.
P. Kellman et al.,
NMRB 2006 19: 352.



Example: Full Reference Scan

- Disable all phase encoding pulses
- All echoes should be equivalent
 - There is T_2^* decay
- Fourier transform in the read direction
 - Shifts \rightarrow phase ramps
- Detect phase difference
- Correct (either phase or delay correction)



Geometric Distortions


- EPI has two readout directions
 - Fast: which is typically considered as readout
 - Slow: which is also named “phase encoding direction”
- Very low readout bandwidth in the PE direction
 - High sensitivity to B0 inhomogeneities ...
 - ... and off-resonance effects such as fat-water shift
- Fat saturation
- Good shimming
- Readout time as short as possible ($BW \sim 1/TA$)




B_0 -induced geometric distortions

- Are shifts in phase encoding direction
- Depend on the phase encoding orientation
- Can be inverted by inverting the PE direction


same
distortion
magnitude!



Gradient Echo




EPI, PE=A-P



EPI, PE=R-L


Reducing distortions with bandwidth




- Maximize PE-bandwidth / minimize echo spacing (ES)

Visualize Distortions in a Phantom


- Phantoms with geometric structures are expensive
 - Structures are likely to induce own distortions
- Use saturation RF pulses to “burn in” lines
 - Multi-slice excitation pulse is a good tool



original 60x60 EPI




mb-sat 60x60 EPI



mb-sat 60x60 EPI pe-inv

Another Reason to Hurry: Resolution Loss



- Shorten full readout time: echo-spacing, BW, parallel imaging

Back to the Basic EPI Readout in Pulseseq

- Plot the k-space trajectory
 - Looks good!? Why worry?
- Let's run `seq.testReport()`

Number of blocks: 121

...

TE: 0.042050 s

TR: 0.084850 s

...

Max gradient:


308e3 280e3 1357e3 Hz/m == 7.25 6.59 31.89 mT/m

Max slew rate:

4.545e9 5.050e9 5.431e9 Hz/m/s == 106.76 118.62 127.57 T/m/s

...

- Let's try to make it faster...



Basic EPI: shorten readout?


ro_duration: 1200us -> 480us

TE: 42.1 ms -> 24.3 ms

ro_duration: 480us -> 240us



TE: 24.3 ms -> 23.2 ms

- Useless “detours” in k-space grow
- Timing does not improve much
- Sampling efficiency suffers




EPI with Ramp Sampling


- Dramatically shorten T_{ESP}
- Increase acquisition efficiency
- Additional steps in the image recon
 - Gridding in the read direction



Basic EPI vs EPI with Ramp Sampling




basic 96x96 EPI



96x96 EPI with ramp sampling

Protocol Optimization and Further Options

- Fat suppression
- Ghosting reduction
- Distortion minimization
 - Distortion correction
- Readout acceleration
 - Multi-shot
 - Parallel imaging
- Overall acceleration – simultaneous multi-slice (SMS)
- Avoid PNS & acoustic resonances




Chemical Shift Artifact & Fat Suppression

- Water-fat separation: ~3.4 ppm
 - At 3T $\sigma \approx 425\text{Hz}$
- Low bandwidth per pixel in PE direction

$$\Delta y = \omega_0 \sigma / BW_{pe};$$
$$425 / 10.4 \approx 40 \text{ pixels}$$



(~30% of FOV)

- Fat suppression
 - Fat saturation, water excitation
- In Pulseseq examples: off-resonance Gaussian pulse before every excitation






Fat Saturation in Pulseseq

- Gaussian pulse at -3.45ppm followed by a spoiler
 - Implementation by some major manufacturer
 - Other vendors use different RF pulses (e.g. minimal phase SLR)




Ghosting: practical advice #1


- For RO use pure gradient
 - Preferably X or Z axis (PNS)
 - Avoid double-oblique slices
 - No in-plane rotation




normal
scaling




overscaled




0°



10°





20°




30°

Ghosting: practical advice #2


- Optimize readout bandwidth / echo spacing
- Do “step left / step right” check
- Avoid mechanical resonances of the gradient system
 - Often the quietest echo spacing setting is the best




1653 Hz/pix




1717 Hz/pix




1786 Hz/pix




1860 Hz/pix



1941 Hz/pix




2029 Hz/pix




2126 Hz/pix

More Troubles to Avoid





mechanical resonances of the gradient system





peripheral nerve stimulation (PNS)

Summary & outlook

- EPI is an established technique
- Good image quality if operated properly
- Numerous new options: SMS, randomized sampling, non-linear image reconstruction
- Need to recognize, understand and avoid artifacts
- Do a pilot run for each specific application
 - Check whether PE inversion gives more favorable distortions



Outlook: towards cutting-edge fMRI



Multi-Vendor SMS-EPI protocol for fMRI applications (ABCD)

<https://harmonizedmri.github.io/>

Acknowledgements:

Berkin Bilgic

Frank Zijlstra

Jon-Fredrik Nielsen

Moritz Zaiss

Qiang Liu

Sebastian Littin

Borjan Gagoski

Imam Shaik

Juergen Hennig

Naveen Murthy

Qingping Chen

Will Grissom

Douglas Noll

Jeff Fessler

Mojtaba Shafiekhani

Niklas Wehkamp

Scott Peltier

Yogesh Rathi

THANK YOU FOR YOUR ATTENTION!



Supported by NIH U24-NS120056 (Nielsen, Zaitsev)
and R01-EB032378 (Rathi, Bilgic)

